

PRODUCT INFORMATION

UNDERWATER TREADMILL



PRODUCT DESCRIPTION



The efficiency of hydrotherapy is widely recognised by medical professionals, and at Variopool too, we are familiar with this popular form of aquatherapy. The use of movable pool floors and bulkheads is a must to guarantee the best results, as they allow the various target groups to make optimal use of the therapy pool facilities, and they ensure the water temperature can be efficiently regulated.

The movable floor can be equipped with movable stairs, a lift for people with disabilities or a ramp to allow rehabilitation patients to access the pool unaided. Underwater treadmills can also be installed, so that patients do not need to turn back once they reach the wall, boosting the efficiency of their workout.

TECHNICAL SPECIFICATIONS AND FEATURES

Material - walking area	Chlorine resistance belt
Speed range (24V DC Drive)	0.8 - 8kph (0.5 - 5mph)
Handrails	Stainless Steel 316 side handrails
Dimensions	1734 x 1184 cm
Training area	L 1405mm (55") x W 855mm (34")
Weight	150 kg
Max distance from main control panel to treadmill	10m
Max distance from main control panel to operator panel	15m
Control system	Touchscreen panel with wireless remote control

SPECIFICATIONS



MORE INFORMATION REQUIRED